




CCA Agar Chromogenic ISO 9308-1:2014

A major revision of the ISO 9308-1 standard has come into force in late 2014. This states that the TTC Agar (Chapman or Tergitol-7 Agar) is replaced by the CCA Chromogenic Agar as a culture medium for the enumeration of *Escherichia coli* and **coliform** bacteria after membrane filtration stage.

The CCA is based on **enzymatic reactions** that give **colour to the colonies** of target organisms for simultaneous detection of coliforms and *E. coli*.

This medium is suitable for **samples with low microbial load** as drinking water, swimming pools, disinfected waters and water treatment plants at the end of the treatment.



Main advantages

- Very good recovery.
- Ideal for detection and collection of *E. coli* and coliforms in water with low contamination.
- Good colour contrast that facilitates the interpretation.
- Available in dehydrated media and 55 mm plates.


Interpretation

E. coli → Colonies from blue to violet


Coliforms → Pink salmon to red

Other bacteria (primarily Gram-negative) →


Colourless (except those with a weak glucuronidase activity but not galactosidase, producing light blue or turquoise colonies).



E.coli



Coliforms



Other enterobacteria
Salmonella



IP-012EN;201702

Technique

1. The water sample is filtered through a membrane filter of 0.45 µm pore diameter, validated according to the ISO Standard 7704:1985 (*).
2. The membrane is then placed on the surface of the CCA medium avoiding entrapment of air bubbles between the membrane and agar surface.
3. The petri dish with the membrane is incubated for 18-24 hours at 36 ± 2°C. If in 18 h there is growth of red or colourless colonies, extend the incubation until 24 h to include late reactions of β-galactosidase or β-glucuronidase.

(*) When the Chromogenic Agar for Coliform is used with the membrane filter method, the colour and growth of the colonies can be modified by the characteristics of the membrane filter. It is advisable to perform a validation of the membrane filter type used.

Composition (g/l)

Enzymatic digest of casein	1.00
Yeast extract	2.00
Sodium chloride	5.00
Monosodium phosphate	2.20
Disodium phosphate	2.70
Tryptophan	1.00
Sodium pyruvate	1.00
Tergitol®7	0.15
Sorbitol	1.00
6-Chloro-3-indoxyl-β-D-galactopyranoside	0.20
5-Bromo-4-chloro-3-indoxyl-β-D-glucuronide	0.10
IPTG	0.10
Agar	13.00
pH: 6.8 ±0.2	

Results

Count β-galactosidase positive colonies and β-glucuronidase negative colonies (all colonies coloured from salmon-rose to red) as Coliform bacteria different from *E. coli*.

Count β-galactosidase positive colonies and β-glucuronidase positive colonies (all colonies coloured from deep blue to violet) as *E. coli*.

Total Coliform count is obtained by the addition of the salmon-rose to red colonies plus the deep blue to violet colonies.

Calculate the concentration of Coliform bacteria and *E. coli* in 100 mL from the initial volume of water filtered and the number of characteristic colonies counted on the membrane. The results are expressed as Colony Forming Units per 100 mL (CFU/100 mL).

Confirmation

To confirm the *E. coli* colonies in this medium a small amount of tryptophan is included verifying indole production:

Coat the blue-violet colonies with a drop of Kovacs' Reagent. If the reagent turns a cherry-red colour in a few seconds this confirms the production of indole and hence the presence of *E. coli*.

Code	Description	Package
417153.1210	CCA Coliforms Chromogenic Agar (ISO 9308-1) (Dehydrated Culture Media)	500 g
447153.0922	CCA Coliforms Chromogenic Agar (ISO 9308-1) (Prepared Plate (Ø 55 mm))	30 plates
Also available CCA Agar according to Spanish regulation SCO/778/2009, of 17 March (not ISO) Ideal for analysis of waste water and very contaminated samples (see technical data sheet)		
446910.0922	CCA Coliforms, Chromogenic Agar (Prepared Plate (Ø 55 mm))	30 plates
Auxiliary Reagent		
252908.1608	Kovacs' Reagent	100 ml

IP-012EN

AppliChem GmbH
Ottoweg 4
DE-64291 Darmstadt
Germany
Phone +49 6151 9357 0
Fax +49 6151 9357 11
info.de@itwreagents.com

Nova Chimica Srl
Via G. Galilei, 47
I-20092 Cinisello Balsamo
(Milano) Italy
Phone +39 02 66045392
Fax +39 02 66045394
info.it@itwreagents.com

PanReac Química SLU
C/ Garraf 2, Polígono Pla de la Bruguera
E-08211 Castellar del Vallès
(Barcelona) Spain
Phone +34 937 489 400
Fax +34 937 489 401
info.es@itwreagents.com

www.itwreagents.com