



pH Adjustment in Pharmaceutical Production

When adjusting the pH in a pharmaceutical manufacturing process it is important to use reagents of the appropriate quality grade. Laboratory products are not suitable for this particular application as some important certificates requested by authorities will be missing.

PanReac AppliChem offers the right product for your needs with the necessary documentation.



Main Advantages

- Acids and bases at different concentrations for different applications (see table 1 and 2).
- Pharma grade: specifications according to main pharmacopoeias (USP and Ph. Eur.).
- Documents and certificates that regulatory bodies require: BSE / TSE, GMO, residual solvents, allergens, trace metal residues, etc.
- Different pack sizes are available to adjust your particular needs.

Table 1: commonly used acids*

Substance	Concentration	Code
Hydrochloric Acid	37 %	141020
	5 M	192109
	1 M	AL1021
Acetic Acid	Glacial	141008
ortho-Phosphoric Acid	85 %	141032
	10 %	147143

*Selection of available substances

Table 2: commonly used bases*

Substance	Concentration	Code
Sodium Hydroxide	Pellets	631687
	6 M	AL6406
	5 M	192159
	1 M	192415
Ammonia	30 %	141130
	25 %	141129

We can develop different concentrations under request. Please feel free to request further information.



IP-018EN;201910

IP-018EN

AppliChem GmbH
Ottoweg 4
DE-64291 Darmstadt
Germany
Phone +49 6151 9357 0
Fax +49 6151 9357 11
info.de@itwreagents.com

Nova Chimica Srl
Via G. Galilei, 47
I-20092 Cinisello Balsamo
(Milano) Italy
Phone +39 02 66045392
Fax +39 02 66045394
info.it@itwreagents.com

PanReac Química SLU
C/ Garraf 2, Polígono Pla de la Bruguera
E-08211 Castellar del Vallès
(Barcelona) Spain
Phone +34 937 489 400
Fax +34 937 489 401
info.es@itwreagents.com



www.itwreagents.com